

FEBRUARY ONLINE SAVERS

VALID ONLINE ONLY UNTIL 29.02.12



GAMACIDE3

- GamaCide3 can be applied to both hard and soft surfaces including: plastic, ceramic, stainless steel, chrome, acrylic, leather, vinyl, glass, Plexiglas®, and baked-on surfaces.
- Non-flammable. No special shipping or storage is required.
- GamaCide3 is non-sticky & leaves no hazardous chemical residue when used as directed. Dye-free.
- GamaCide3 is compatible with most impression materials and dentures.
- GamaCide3 is biodegradable.
- GamaCide3 can safely be disposed via standard waste disposal methods.

\$20.00 EA



SHIELD FORCE PLUS

Tokuyama® Shield Force Plus is a one-component, single application, light-cured protective sealant designed for the treatment of hypersensitive dentin. Tokuyama® Shield Force Plus effectively infiltrates and seals exposed dentinal tubules to inhibit the flow of fluid through the tubules which is the cause of sensitivity. At the same time the durable coat of cured Tokuyama® Shield Force Plus helps prevent loss of dentin due to abrasion and/or erosion of exposed cervical dentin.

- Quick and effective treatment of hypersensitive dentin

- Simple procedure
- Durable coating on exposed dentinal tubules
- Superior sealing of prepared tooth
- Treatment of hypersensitive dentin
- Reduction of abrasion and erosion of exposed cervical dentin
- Alleviation and/or prevention of tooth sensitivity after tooth preparation for direct and indirect restorations

\$99.00 REFILL



INJEX

**LOCAL ANAESTHESIA WITHOUT NEEDLE.
INJEX: THE SOFT WAY TO INJECT.**

The fear of needle injections has been widespread since the development of the conventional syringe in 1853. Experienced dentists are well aware of the different problems associated with the conventional needle syringe: Delay until onset of anaesthesia, painful injections, feelings of fear and anxiety exhibited by patients,

tissue damage, and unpleasant numbness of soft tissue, risk of injury with restless or fearful patients...

The INJEX System uses an injection ampoule with a micro orifice of only $\varnothing 0.18$ mm through which the anaesthetic is administered under dosed pressure to the submucosa – virtually painless and exactly where it is needed.

\$350.00 EA



DIATECH DIAMOND KIT FOR CROWN & BRIDGE

DIATECH Multilayer diamonds are multiple layered in an exclusive manufacturing process using the most advanced technology, while the diamond instruments you used in the past probably had just one layer of diamond. Multiple layered diamonds are provided with more diamond on different levels and in various layers. New diamond grits with sharp edges will, therefore, take the place of broken off grains or grains separated from the plating during normal use. DIATECH Multilayer diamonds have an extended life and cut faster.

\$199.00 KIT

